

# Napa™

PERFECT CHOICE Carpet Cushion



## Luxurious Cushion Comfort

Like walking on a cloud

**PET STOP**  
TECHNOLOGY™

Block-Out® Barrier  
Bottom Moisture Barrier  
Fresh Dimension® Antimicrobial

## Perfect Choice for Plush and Textured Carpets

### Comfort and Durability

- DuraPlush® foam technology adds superior cushioning and a rich fullness to your carpet
- Condensed, pure urethane foam ensures durability and long-lasting beauty of carpet

### Pet Stop Technology™

- Fresh Dimension® antimicrobial inhibits the growth of mold, mildew, and odor-causing bacteria
- Block-Out® moisture barrier repels spills
- Bottom moisture barrier provides additional protection

### Healthy Home

- Hypoallergenic
- CRI Green Label Plus-certified for indoor air quality

### Fundamentally Green

- 100% recyclable at the end of its useful life
- Contributes toward achieving LEED® building credits

### Peace of Mind

- Insulates home for better energy efficiency
- Reduces noise for a quieter, more enjoyable home
- Life-of-the-Home warranty with additional StairGuard™ protection


Learn more at [LPFlooringProducts.com](http://LPFlooringProducts.com)



*Leggett & Platt*  
FLOORING PRODUCTS

# Napa™

PERFECT CHOICE Carpet Cushion



## SPECIFICATIONS

|                      |  |
|----------------------|--|
| Product code:        | BZ0288   |
| R-value:             | 2.22   |
| Thickness:           | 0.500"   |
| FHA class:           | 2  |
| Roll size:           | 6' x 45' (30 yd <sup>2</sup> )   |
| Cushion color:       | Blue   |
| Pet Stop Technology: | - Fresh Dimension® antimicrobial<br>- Block-Out® top moisture barrier<br>- Bottom moisture barrier |
| Certification:       | - FHA/HUD/VA (UM72a per ASTM D3574)<br>- CRI Green Label Plus                                      |



## RECOMMENDED INSTALLATION INSTRUCTIONS

1. Always inspect the floor prior to installation to ensure it's clean and free of debris. **Inspect and test for moisture before installing any dual-moisture-barrier product.** Reference the CRI Carpet Installation Standard.
2. When laying cushion over subfloor where old glue is exposed, first sprinkle talcum, baby powder, or baking soda. This will help prevent the cushion from partially sticking to the floor, potentially causing noise.
3. If using a pad adhesive, apply a small amount in the corners of the room to help hold the cushion in place during installation. Staples may also be used around the perimeter of the room.
4. Use a flexible, moisture-resistant-type tape (no masking or other paper-based tape) at cushion seams and for covering tears, punctures, and staples.

## LIFE-OF-THE-HOME LIMITED WARRANTY

Please visit [LPFlooringProducts.com/napa](https://LPFlooringProducts.com/napa) for detailed warranty information.



*Leggett & Platt*  
FLOORING PRODUCTS